


S@FE SP Special Cabinets for Robot



S@fe SP cabinets have been designed in order to provide a controlled environment in which place liquid handling robots. The cabinets provide the same level of containment of a Class II cabinet as defined by the EN12469:2000 norm ($Apf > 10^5$).

With up to 2.4m of internal width available, S@fe SP cabinets provide enough room to comfortably fit liquid handling robots.

In order to allow easy installation of the cabinet and the robot, S@fe SP are designed to be completely disassembled before shipping and reassembled in place at customer site!

