

# CELLCARE

CO<sub>2</sub> INCUBATOR  
EXCELLENCE FOR YOUR LAB



## PRECISION & COMFORT

- The Heart - ensure stable humidity and CO<sub>2</sub> levels and accurate readings with a multi-parameter sensor
- O<sub>2</sub> control and regulation system available on request
- The Brain - keep your cells safe with an advanced electronic control system
- Real-time data storage with a user-friendly touch screen interface with USB and ethernet connectivity
- No internal spare parts to handle (e.g. UV lamps or HEPA filter)



## CLEANABILITY & CONTAMINATION CONTROL

- 180° C (356° F) decontamination
- Ensure optimal cell growth thanks to a seamless deep drawn chamber with fanless design
- Removable and autoclave-safe shelving system and water tray
- Smooth external surfaces with painted steel

### Dimensions

53L  
188L


## PERFORMANCE & SPEED

- Proprietary advanced direct heating system for temperature uniformity
- CO<sub>2</sub> and temperature recovery <5 minutes
- Protect cells from vibration and turbulence with a seamless design



## SAFETY

- No need to remove sensor during high temperature cycle
- Anti-vibration soft door with available automatic locking system



# DIMENSIONS

| MODEL                    | 188                                | 50                                  |
|--------------------------|------------------------------------|-------------------------------------|
| External dimensions (in) | 27.79(W) x 30.31(D)<br>x 43.9 (H)  | 24.01(W) x 24.41(D)<br>x 37.07(H)   |
| Internal dimensions (in) | 20.87(W) x 19.68(D)<br>x 27.16 (H) | 15.75(W) x 13.78 (D)<br>x 15.75 (H) |
| Volume (L)               | 188.6                              | 53                                  |

FIND OUT MORE



YOUR SAFETY  
IS OUR COMMITMENT

**TECNILABO**  
TECNIPLAST GROUP

TECNILABO HEADQUARTERS

TECNIPLAST S.P.A.

VIA I MAGGIO, 6

21020 BUGUGGIATE (VA)

tel +39 0332 809711 - fax +39 0332 458315

info@tecnilabo.com

SISTER COMPANIES

CANADA | Toll Free: +1 855.347.8718 • [www.tecniplastcanada.com](http://www.tecniplastcanada.com) • [info@tecniplastcanada.com](mailto:info@tecniplastcanada.com)

CHINA | Tel. + 86 (0)21 50810920 • [www.tecniplast.cn](http://www.tecniplast.cn) • [tecniplastchina@tecniplast.it](mailto:tecniplastchina@tecniplast.it)

FRANCE | Tél. +33 (0)4 72 52 94 41 • [www.tecniplast.fr](http://www.tecniplast.fr) • [info@tecniplast.fr](mailto:info@tecniplast.fr)

SWITZERLAND | Tel +41 (0)81 750 31 40 • [www.tecniplast.ch](http://www.tecniplast.ch) • [info@tecniplast.ch](mailto:info@tecniplast.ch)

UNITED KINGDOM | Tel. + 44 (0)345 0504556 • [www.tecniplastuk.com](http://www.tecniplastuk.com) • [info@tecniplastuk.com](mailto:info@tecniplastuk.com)

USA | Toll Free: +1 877.669.2243 • [www.tecniplastusa.com](http://www.tecniplastusa.com) • [info@tecniplastusa.com](mailto:info@tecniplastusa.com)